



Mestrelab Research

chemistry software solutions

Mnova qNMR 2.0 – the tool for all skill levels

Mike Bernstein, VP R&D

28 March 2018



mike@mestrelab.com



+44 7910 371511


RESOURCES

[LATEST](#) [ARTICLES & NEWS](#) [PUBLICATIONS](#) [RESOURCES](#) [CHANGELOG](#) [VIDEOS & WEBINARS](#)

YOU ARE AT: [Home](#) » [Search Results for "qNMR" Query](#)

SEARCH RESULTS: qNMR (50)

VIDEO [Processing and analysis in Mnova 12. focus o...](#)




BY ENRIQUE SÁNCHEZ — 14 FEBRUARY, 2018

Processing and Analysis in Mnova 12. Focus on requirements for qNMR

In this video Dr. Mike Bernstein focuses on requirements for qNMR analysis using Mnova 12....

MNOVA QNMR



BY ENRIQUE SÁNCHEZ — 14 FEBRUARY, 2018


Processing and Analysis in Mnova 12. Focus on requirements for qNMR (pdf)

Title: Processing and Analysis in Mnova 12. Focus on requirements for qNMR Author: Mike Bernstein...

WEBINAR [R - 2018 Webinar Series Requirements for NMR Internal standards](#)



MNOVA QNMR



SEARCH

UPDATES

2018

MAR 23

What's new in Mnova 12.0.2

JAN 09

What's new in Mnova 12.0.1

2017

OCT 11

What's new in SMA 2.0

NEWSLETTER SUBSCRIPTIONS



You **Tube** → Mestrelab channel



<http://www.nmrvalidation.org/>

qnmr.org



Agenda

qNMR basics

Results expectations

Do's and don't of processing

Integration options

qNMR Purity plug-in - *DEMO*

SMA to perform purity – *DEMO*


Automated Purity – *DEMO*

Questions



The fundamentals are well known

J Shoolery, "A Basic Guide to NMR", 1972



Internal reference method



Test sample

known compound(s)' weight(s) known
→ purity to be determined


Quantitation Reference Material



known purity

- no overlapping signals in a given matrix
- precisely weighable (non-hygroscopic / non-volatile)
- chemically inert with respect to the solvents used and the matrix analyzed






Scheme 1.1. Interrelations between the concepts true value, measured value, error and uncertainty.

Complete uncertainty analysis:

- Trueness
- Precision
- Selectivity
- Limit of detection
- Limit of quantification
- Linear range
- Working range
- Robustness
- Measurement uncertainty



How accurate and precise do you need your qNMR result?



“High performance qNMR” ($u < 0.5$)



Accurate purity ($u < 2.0$)



Standard purity ($u < 5.0$)

“Quantitative” NMR spectrum

Full peak integrals

Good peak-shape

Signal-noise-ratio (SNR)

Time between excitation pulses (“relaxation delay”)

Temperature control of the sample

Shimming

Total data acquisition time (number of scans)

Spectrum processing

Apodisation/Line-broadening

improves peak shape at base (“wiggles”)

Data points for FT (zero filling)


for effective integration

Phase correction

for effective integration

Baseline correction

for effective integration




“Processing and analysis in Mnova 12. focus on requirements for qNMR”


<https://www.youtube.com/watch?v=cctl2I2XjwQ>

“Regions” autophase

Automatic phasing can provide consistency, and often outperform manual phasing



Polynomial baseline correction



Processing templates

Use this for *consistency*

Processing Template

f1 More Processing Analysis

Time Domain

☐ Truncate

☐ Drift Correction

☐ FID Shift

☐ Frequency Shift

☒ Apodization
Exponential: 0.25 Hz

Zero Filling and LP
Spectrum Size: 65536

Fourier Transform
Protocol: None
Swap Halves: on
Mirror Image: on

Frequency Domain

☒ Phase Correction
Method: Regions Analysis

☒ Baseline Correction
Method: Bernstein Polynomial Fit
Polynomial Order: 3

☐ Smoothing

☐ Reverse

☐ Reference

☐ Absolute Reference
Proton Reference: 400
Factor: 100.000000

☐ Cuts
Number of Cuts: 0

Apply OK Cancel

Integration method – IMPORTANT!



GSD – good when overlap or close multiplets, but precision/accuracy is worse



Sum – this is the “standard” method



qGSD (\geq Ver 12.0) – the accuracy of Sum, and handles close or overlapping peaks




Edited sum – improved Sum

Mnova qNMR uses:

Standard integration → for the reference peak(s)

Multiplets integration → for the solute multiplets

Good



Desirable but not totally necessary

Chemical structure of analyte

Better nuclides determination (AutoAssign)

Better peak “types” (impurity, etc.)

Sample metadata in the title

Good laboratory practice

Transcriptional errors ↓

Metadata co-exists with the data – for ever



00110779-990-b_2011096114309

SW:4.518

MW:410.471

RW:1.906

CR:DMP

How does Mnova do Purity qNMR (1-3: Purity plug-in)?

Reference materials' essential data are stored (and can be shared)

Sample spectrum processing templates

“One off” new sample analysis **(1)**, or fixed SOP “experiment” **(2)**

Replicates and repeats **(3)**

Mnova will select the best multiplets for quantitation

- How many to select and average
- Selection “rules”

Multiplet integration limits are automatically set

Nuclides *per* multiplet are determined (→ “Normalized” integrals)

Automatic selections can be manually changed

- Nuclides
- Integral range

How does Mnova do Purity qNMR (4: SMA)?

A flexible, quantitation tool that you can tailor to your needs

1D and 2D regions integration

Automatic equation generation, including qNMR Purity

How does Mnova do *Automatic Purity* qNMR (5: Batch qNMR)?

Based on conventional Mnova qNMR Purity plugin

Mgears* automation-

Input – batch, listener, DB, command line interface

Processing file

Analysis, e.g. Purity

Output format

Output file(s) to be written

*Mnova Ver 12.0.2 and higher

Conclusions

Mnova has strong, “basic” processing and analysis capabilities

Purity analysis can be easily and reliably performed

- single sample

- repeats or replicate

- prescribed analysis

More complex quantitation is performed using SMA

Full-automated operation can be performed

- e.g., JEOL resonance “qNMR Seamless”