

<b><i>ERROR NAME</i></b>	<b><i>CODE</i></b>	<b><i>DESCRIPTION</i></b>
E_LIC_OK	0x00000000	There is no error (This code will never be displayed in the error message dialog box, as it means that everything went ok)
E_LIC_FATAL	0x80000001	Fatal error (This should not be a frequent error).
E_LIC_TIMEFAKED	0x80000002	It was detected that the user has set back the system clock
E_LIC_EXPIRED	0x80000003	The license has expired (Expiry Date exceeded)
E_LIC_INVALID	0x80000006	Invalid license (most likely due to an error in generating the license)
E_LIC_UNKNOWN	0x80000007	Default value. It is displayed when a license is not in use (has a blue circle with white question mark)
E_LIC_NOTAVAIL	0x80000008	There are no more available licenses. This error is displayed when a client with a campus or floating license tries to obtain a license from the Server and there are no more left.
E_LIC_KEYNOTFOUND	0x80000009	This message should not appear, as it indicates that the public key could not be loaded, necessary to verify the license. If this occurs, then there is a strange error in the encryption module of the OS.
E_LIC_NOTFOUND	0x8000000A	A license for the product was not found (This message normally appears when there is no license)
E_LIC_INVALID_SIGN	0x8000000B	Invalid license signature (Indicates that the user tried to modify the license file, or when sending the license by email, the license file was modified in some manner)
E_LIC_INVALID_HOST	0x8000000C	Invalid Host Id (Indicates that the Host id does not correspond to the one included in the license)
E_LIC_INVALID_VERSION	0x8000000D	Invalid product version (E.g. The license is for NMR v1.0 but the user installed NMR v2.0 and has an expired U&S)
E_LIC_TEMP_RIGHTS	0x8000000E	The program could not access the <i>cache</i> folder because most likely does not have the right permissions to do so
E_LIC_SRV_CONNECT	0x8000000F	This is a <i>floating</i> or <i>campus</i> license which could not connect to the server. There are several reasons: <ul style="list-style-type: none"> <li>• The string &lt;server&gt; was not configured in the license file</li> <li>• The firewall does not allow the connection</li> </ul>
E_LIC_SRV_PROTO	0x80000010	These errors are normally due to temporary

<b><i>ERROR NAME</i></b>	<b><i>CODE</i></b>	<b><i>DESCRIPTION</i></b>
E_LIC_SRV_BADREQUEST	0x80000011	failures in the connection to the server, and normally, trying again resolves the problem. This error may also mean that an incompatible version of the server is being used
E_LIC_SRV_BADRESPONSE	0x80000012	Due to temporary error in the connection to the server -or- the host ID of the server does not
E_LIC_SRV_BADSIGN	0x80000013	Infrequent error, due to either a strange failure in the network, or someone tried to replace the license server
E_LIC_SRV_BADDATE	0x80000014	The date of the license server and the client license differs in more than 2 days. One of the dates is incorrect
E_LIC_SRV_SESSIONNOTFOUND	0x80000015	This error is corrected upon restarting the program. It means that the server could not find the session for the client
E_LIC_SRV_DENY	0x80000016	The administrator of the server explicitly banned this client from obtaining a license (the license server administrator must be contacted directly)